





Preference Menu

The Preferences window allows for customization of the Sonata® System. When an item on the list is chosen, it will be displayed and available for selection from the Ultrasound Display.



Using the Preferences menu, you can:

- Add or remove optional functions from the Ultrasound Display.
- Set System Date and Time.
- Choose Language for the User Interface.
- Select the preferred Controls:
 -  : Display Ultrasound Controls for image optimization
 -  : Hide SMARTTM Guide or temporarily hide graphics
 -  : Show Scale or approximation of size
- Manage User Profiles under the System tab.



For more information, refer to the ***Sonata Operator's Manual*** or the ***Quick Reference Guide***.